

THE BITE SPLINT WITH THERMO-EFFECT



**excited
patients**

**pleased
dentists**

**relaxed
technicians**

PREMIOtemp CLEAR FLEX

Nothing presses, nothing straints

**the milled PREMIOtemp CLEAR FLEX bite
splint is just super comfortable to wear.**

PATIENT INFORMATION



CLEAR FLEX

MAXIMUM BREAKAGE RESISTANCE WITH THERMO-EFFECT

REFORMO



Maximum breakage resistance at body temperature through thermoplastic flexibility

TOP FEATURES

YOUR BENEFITS



- Tension-free wearing comfort and precise fitting due to the thermo-effect



- Very durable – maximum breakage resistance through memory effect – more resistant to breakage than most conventional PMMA plastics
- High transparency for almost invisible wearing
- Very low plaque affinity

Memory-Effect: 100% reset at body temperature without loss of surface hardness

CLEAR FLEX

CLEAR FLEX

**The milled PREMIOtemp
CLEAR FLEX bite splint from primotec**

Care instructions:

- The PREMIOtemp CLEAR FLEX splint should be cleaned under running water with toothpaste after each use
- Store your bite splint in a dry and dark place if not worn
- Do not use denture cleaners / alcohol-based cleaners
- Place the bite splint in warm water at 40-50 °C for 1-2 minutes before insertion
- Protect the splint from extreme heat and coldness



Joachim Mosch e.K.
Tannenwaldallee 4
D-61348 Bad Homburg
Tel. +49 (0) 61 72 - 99 77 0 - 0
Fax +49 (0) 61 72 - 99 77 0 - 99
primotec@primogroup.de
www.primogroup.de